



Cleaning procedure for ExiGo flow sensors

TABLE OF CONTENTS

1.	General handling	3
2.	Examples & Cleaning Procedures	4
3.	Cleaning Methods that are not Recommended	7

1. GENERAL HANDLING

Do not allow the sensor to dry with media in the capillary tube without flushing clean first. Also try to avoid letting the filled sensor sit for extended periods (depending on your liquid).

Before storing the sensor, always drain of fluid, flush with cleaning agent, blow out, and dry the capillary. Sensor plugs should be installed for storage. The cleaning agent (detergent, solvent, etc.) should be chosen for its effectiveness in removing the liquid media, and compatibility with borosilicate glass.

In general we recommend finishing any cleaning with an acetone and/or IPA flush before letting the sensor dry out and closing it for storage. The advantage of using clean IPA is that it evaporates without residue (in contrast to acetone, for example).

2. EXAMPLES & CLEANING PROCEDURES

- **Working with Multiple Liquids:**

Switching between multiple liquids can leave transient deposits in the form of liquid layers inside the glass capillary. This is especially common for insoluble liquids, but can happen even with miscible liquid combinations. For example, when IPA is followed by water in a sensor without drying in between, large offsets can be observed for hours after switching to water.

If possible, dedicate a separate sensor for each different liquid to be measured. If not possible, use caution when switching media and clean properly.

- **Working with Water:**

When working with water it is recommended to **not** let the sensor dry out. All salts and minerals in the water will deposit on the glass and are difficult to remove. Although salt solutions are particularly prone to problems, even clean water can still contain enough dissolved minerals to form a deposition layer. Flush with DI water on a regular basis to prevent build-up. If you still encounter problems, occasionally flush the sensor with slightly acidic cleaning agents.

When working with water containing organic materials (sugars, etc.) microorganisms often grow on the walls of the glass capillary and form an organic film that can be difficult to remove. Flush on a regular basis with solvents such as acetone, ethanol, methanol or IPA, or with cleaning detergents to remove organic films.

- **Working with Oils:**

When working with oil it is recommended to **not** let the sensor dry out. Edible oils can become rancid, dry out or simply change their composition over time and leave sticky films on the wall of the capillary. This will typically happen if the sensor is only drained after usage but not properly cleaned. Then the oil film on the capillary wall becomes a constant deposition and will be more difficult to clean than right after emptying the sensor. These films might not be dissolved by oil anymore and cause measurement deviations.

Flush with solvents such as acetone, methanol, ethanol, IPA, etc or cleaning detergents on a regular basis to remove oil films.

- **Working with Silicone Oils**

When working with silicone oil it is recommended to **not** let the sensor dry out. Silicone oils can be cleaned out using acetone or special cleaners. Check with your silicone oil supplier for cleaning agents compatible with glass surfaces.

- **Working with Paints or Glues**

When working with paints or glues it is critical to **not** let the sensor dry out. Often, depositions of paints and glues cannot be removed anymore after they have dried. Flush the sensor with cleaning agents recommended by your paint or glue manufacturer that are compatible with glass. Ensure that you have found a good cleaning procedure before performing the first tests, and always clean shortly after emptying the sensor.

- **Working with Alcohols or Solvents**

Unlike most other fluids, alcohols and solvents are not critical and a short flush of acetone followed by IPA (or only IPA) is sufficient to clean the capillary walls.

- **Other Liquids or Applications**

If uncertain about your application and how to clean the flow sensor, please contact Cellix for additional support

3. CLEANING METHODS THAT ARE NOT RECOMMENDED

In general, **any cleaning by mechanical means should be avoided.**

Never enter the sensor's flow path with sharp objects that could scratch the glass surface.

Furthermore, no abrasives or liquids containing solids that can grind the surface clean should be used. Anything that affects the glass wall will cause deviations in the measurement performance or permanently damage the sensor.

Strong acids and bases should also not be used to clean the sensor. Acids can sometimes be used in low concentration and at low temperatures. Before using the acid check how compatible it is with borosilicate 3.3 glass (Pyrex® or Duran®).

For more information on the ExiGo Pump or any other Cellix product or service, please call:

Republic of Ireland: +353 (0) 1 4500 155

This publication is not intended to form any contract for computer software with Cellix Limited or any of its subsidiaries. If you require a particular application or feature, the suitability of the program should be verified.

Every effort has been made to ensure the contents of this and all Cellix publications are accurate and correct. This publication is for information purposes only and is not intended to form the basis of any contract for the sale or purchase of the product being described. No liabilities for errors will be accepted.



CELLIX Limited

Unit 1, Longmile Business Park,

Longmile Road,

Dublin 12,

Ireland

Tel: +353-1-4500-155

Fax: +353-1-4500-158

Email: info@cellixltd.com

Web: www.cellixltd.com

© Cellix Ltd. 2005 – 2015